


Materials Compatibility of the Chimney Camera Intrusions -Issues and Plans -

K. Heeger, LBNL

Chimney Monitoring Cameras



NOTE: CAMERA SITS IN THE SAME AIR AS THIS PERSON BREATHEES.



Purpose:


To visually monitor the deployment and retraction of the 4pi calibration system.

10 m

Note: only spheres, cylinders, and deck level are accurately scaled.
Positions of straps, PMTs, light shield, camera, and upper balloon shape are schematic only.

Materials Issues


NOTE: CAMERA SITS IN
THE SAME AIR AS THIS
PERSON BREATHEES.



To be immersed in LS, or not?
What are long-term effects?


Standard MDC window
(covlar glass-to-stainless weld)
Custom MDC window
(gold glass-to-stainless weld)

Or, custom design?




Step 1

Deploy camera intrusion into gas volume only and try it. If we need to deploy into LS we can always extend tube.




Note: neck PMTs provide precedence for intrusion into gas volume.



Step 2 - Option 1

Deploy camera intrusion into LS volume with MDC window and welded flange.


-> requires long-term LS test for MDC window.



Step 2 - Option 2

Deploy camera intrusion into LS volume with custom lucite enclosure.
(Kishimoto's idea)

-> requires machining of custom lucite enclosure



The Next Steps

- Have already procured standard MDC window, will start soak test soon. (LBNL+Mozumi)
- Will inquire about availability of custom MDC window with gold weld. (LBNL)
- Complete design of camera intrusions and provide drawings. (LBNL)
- Select cameras and develop readout system. (KSU)
- Perform illumination tests. (LSU)

Final Note: Chimney monitoring cameras will be very useful but should not hold up installation of the 4pi system.